



Portable Pressure Vessel & Vessel Systems



Reagent Delivery System - 30L with 12" Diameter Holding Vessel on Support Stand with sensors on Main Cart-Shown Attached and Detached.



Portable Processing Systems complete with motors, support stands with casters, with mounted load-cells - consists of one 60 liter vessel 12" diameter with removeable head and 140 liter vessel with removeable head.



Portable Pressure Vessel Compounding System, made of stainless steel with a 30 Gallon (100 Liter) vessel with offset homogenizing assembly and variable speed stirrer assembly. Complete with a flanged bottom, load cells, tachometer, the unit sits on a stainless steel cart complete with casters, handle and brakes.

Compounding vessels to consist of 25 liter and 50 liter postable suspension vessels with mixing systems. Mounted on carts with pneumatic device to raise and lower removable heads.

GO TO PPV IN ACTION:

